

Choose your **C Line** *touch* Plan

Gold Plan

£3190 + VAT

C Line Touch Machine + 10 vials (10 treatments) + Full Training including marketing assistance and package (on SEO, Social Media)

Platinum Plan

£3490 + VAT

C Line Touch Machine + 30 vials (30+ treatments) + 1Need Micro-needling Pen + (30 cartridges) + Full Training including marketing assistance and package (on SEO, Social Media)

1NEED



C Line *touch*



HYA - TOUCH FACIAL



0746 339 39 29



info@revitacare.uk



www.revitacare.uk

LABORATOIRE
REVITACARE
FRANCE

IRB TRADE | Revitacare UK

C Line

touch

HYA - TOUCH FACIAL

A multi step facial incorporating Electroporation and the latest Micro-current technology allowing product penetration of the Highest Quality products made by Laboratories Revitacare a french laboratory founded in 2003, specialised in research, design, manufacturing of medical devices based on hyaluronic acid.

This Machine allows for upwards of 36 different combinations via 3 different modalities to provide your clients with the bespoke tailored facial they would need! Allowing you to treat all skin types and ages.

The C Line Touch Machine has been specifically designed to cater to the molecular weight of Revitacare products providing 100% of the product to penetrate into the skin without using any needles.

Enquire today!



0746 339 39 29



info@revitacare.uk



www.revitacare.uk

C Line

touch

HYA - TOUCH
FACIAL

Price List

Ex VAT

Sterile Cosmetics Range	Single Box	5+ Boxes
Cytocare 715 C LINE (5x5 ml)	£145	£135
Cytocare 640 C LINE (5x4 ml)	£120	£110
HairCare C Line (10x5 ml)	£125	£115
CelluCare C LINE (10x5 ml)	£125	£115
StretchCare C LINE (10x5 ml)	£125	£115



0746 339 39 29



info@revitacare.uk



www.revitacare.uk

HYA - TOUCH FACIAL

C Line *touch*

HYA - TOUCH FACIAL

The **C Line Touch** machine offers a non-invasive technology with an electroporation program developed to optimize product penetration.



La Beauté commence par la qualité de la peau. 173.
"Beauty starts with skin quality."

HYA - TOUCH FACIAL

Non-invasive and painless, the C Line Touch invents the manual beauty care of tomorrow by integrating unique power of Electroporation and micro-currents in a session.- Giving you a unique facial everytime!

Rejuvenate Program

- A combination of currents to boost the cosmetic actives according to the skin's needs.
- Anti-ageing treatments
- Active ingredient Customisable program

Glow Program

- An exclusive light current to improve blood microcirculation and blood flow.
- Dark circles, dull skin Venous dilatation especially the eye area!
- Increase in oxygen levels - for added oxygenation to the skin

Lift program

- A stimulating current for the Electrolift of the skin muscles and the face.
- Stimulation Toning
- Skin firming



Specifically Designed to the molecular weight of Revitacare products to give you 100% product penetration into the cells!

☎ _____
0746 339 39 29

✉ _____
info@revitacare.uk

🌐 _____
www.revitacare.uk



Totally Wireless -Equipped with advanced technology, USER ASSIST®

Bluetooth connectivity

Ideal comfort in session. The C Line Touch works like an exoskeleton on your fingers. It stands out for its precision in sensitive areas thanks to its manual and intuitive use. Combines the agility of your movements with the intuitiveness of your grip to guarantee a pleasant session comfort. Eye contour, lip contour with ease.

Instantly visible results

The C Line Touch offers universal penetration of cosmetic ingredients and immediate results.

Exclusively, the C Line Touch benefits from the expertise of our engineers with the intelligent software USER ASSIST®. Equipped with advanced technology, USER ASSIST® is a real assistant in your daily life and a force of proposal in your sessions. USER ASSIST® analyses the needs and information of the skin in order to create a treatment protocol that is personalised to the expectations of each client. In fact, this software instantly programs all the session parameters (duration, microcurrent, frequency used, etc.) thanks to its numerous reference treatments.

Thus, the collection of this skin data by USER ASSIST® configures the ideal protocol for each client and ensures the best possible results.


So What is Electroporation?

Electroporation is the application of an electric current passing through the skin and cell membrane in order to open pores and create channels through which product can penetrate into the skin cells and to the dermis. The C-Line Touch is an intradermal electroporation device that penetrates the epidermal and dermal layers of the skin allowing for up to 100% penetration of the product! Allowing more efficient delivery of products into the skin cells. Specifically designed to penetrate Revitacare C-Line products. This Technology can be used for the face and body including sensitive areas such as the eyes and lips!

The applicators are designed so that they can reach under the eye and target smaller areas where larger devices cannot.


Electroporation can help:

- Expression lines and wrinkles
- Face firming and toning
- Anti-aging
- Hyper-pigmentation
- Rosacea
- Oily skin
- Skin rejuvenation
- Skin hydration
- Brightening of the skin
- Lifting



Hya-Touch Facial

Designed to include the perfect facial for your clients to Glow- Enquire today!

 [0746 339 39 29](tel:07463393929)

 info@revitacare.uk

 www.revitacare.uk